



Calf Intestinal Alkaline Phosphatase (CIAP)

Cat. No. 18009-019

Size: 1,000 units

Lot No.

Conc.: U/μl

Exp. Date:

Store at -20°C.

Description:

NOTE: DILUTION BUFFER FORMULATION HAS BEEN CHANGED. THIS PRODUCT IS NO LONGER COMPATIBLE WITH 18059-014.

Calf intestinal alkaline phosphatase (CIAP) is a phosphomonoesterase purified from calf intestinal mucosa that hydrolyzes 5'-phosphate groups from DNA, RNA, and nucleotides. CIAP is used to dephosphorylate linearized vector DNA prior to insert ligation and to remove 5'-phosphate groups prior to 5'-end labeling of nucleic acids with T4 polynucleotide kinase.

Components:

18009-019	CIAP	Lot No
Y01371	CIAP 10X Buffer	Lot No
50839	CIAP Dilution Buffer	Lot No.

Unit Definition: One unit of CIAP hydrolyzes 1 mmol of p-nitrophenyl phosphate in 1 min at 37°C.

Dilution Buffer:

25 mM Tris-HCl (pH 7.6)
1 mM MgCl₂
0.1 mM ZnCl₂
50% glycerol (v/v)

Unit Assay Conditions

1 M Diethanolamine (pH 9.8)
0.25 mM MgCl₂
10 mM p-nitrophenyl phosphate

Dephosphorylation Buffer (1X concentration): 50 mM Tris-HCl (pH 8.5), 0.1 mM EDTA.

The Dephosphorylation buffer is supplied as the 10X concentrate and should be diluted, 1:10 (1 part Dephosphorylation buffer + 9 parts of other components = 10 parts final reaction volume.)

Quality Control Assays:

No detectable contaminating activity is observed in endodeoxyribonuclease, exodeoxyribonuclease, and ribonuclease assays. The enclosed buffers were assayed with the enzyme and met quality control specifications.

Doc. Rev. 091499

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

For technical questions about this product, call the Life Technologies TECH-LINE™ (800) 828-6686.

Traditional Protocol:

This protocol dephosphorylates 1 pmol of 5'- DNA termini from purified DNA. DNA dephosphorylated by this method is suitable for cloning or for labeling by T4 polynucleotide kinase using the Forward Reaction:

1. Determine the mass of DNA required for 1 pmol of the type of DNA 5' end.
2. To a 1.5-ml microcentrifuge tube, add 4 μ l of CIAP 10X Buffer and 1 pmol of DNA ends.
3. Add autoclaved, distilled water to 39 μ l.
4. Dilute CIAP in dilution buffer such that 1 μ l contains the amount of enzyme required for the appropriate 5' end (*i.e.*, 1 unit for 5'-recessed and blunt ends and 0.01 units for a 5' overhang).
5. For 5'-recessed and blunt-ended DNA, incubate at 50°C for 60 min. For DNA with a 5' overhang, incubate at 37°C for 30 min.
6. Inactivate/remove the CIAP according to the protocol described below.

Simplified Protocol:

This protocol allows for the dephosphorylation of DNA directly in restriction endonuclease buffer in the presence of the restriction endonuclease. This is a convenient way of preparing DNA for cloning.

1. Restriction endonuclease digest the vector DNA. (**NOTE:** Heat inactivation of the restriction endonuclease and subsequent purification of the vector DNA are not necessary.)
2. Add 1 unit of CIAP to the restriction endonuclease digest.
3. For 5'-recessed and blunt-ended DNA, incubate at 50°C for 5 min. For DNA with a 5' overhang, incubate at 37°C for 5 min.
4. Inactivate/remove the CIAP according to the protocol described below.

Inactivation/Removal of Calf Intestinal Alkaline Phosphatase (3 methods):

1. Heat Inactivation: Note the $MgCl_2$ concentration in the reaction and add EDTA (pH 8.0) to an equal final concentration. Incubate the reaction at 65°C for 15 min.
2. Organic Extraction: Add an equal volume of phenol:chloroform:isoamyl alcohol (25:24:1). Vortex thoroughly and centrifuge at $14,000 \times g$ at room temperature for 5 min. Carefully remove the upper, aqueous phase and transfer it to a fresh microcentrifuge tube. Add 0.1 volume of 3 M sodium acetate. Vortex. Add 2.5 volumes of 100% EtOH. (**NOTE:** Do not substitute NH_4OAc for NaOAc because NH_4 ions inhibit T4 polynucleotide kinase.) Vortex the mixture thoroughly and centrifuge at $14,000 \times g$ at room temperature for 5 min.
3. CONCERT™ Gel Extraction Systems: Following electrophoresis of the dephosphorylated DNA on an agarose gel, use the protocols supplied with the CONCERT Gel Extraction Systems to purify the DNA.